## Separation and Analysis of Three Unknown Solids

Steven Labalme Lab Section 1A05

November 19, 2021

Lab Assignment 1d CHEM 22000

## **Figure**




Figure 1: Isolation and purification process for a selection of unknown solids.

Lab Assignment 1d CHEM 22000

## Lab Report

- 1. What is the identity of unknown A and how did you determine this?
  - Unknown A is resorcinol. The experimental melting point range of 108.7 °C to 109.6 °C falls within the known melting point range of resorcinol, and is disjoint from all other known melting point ranges.
- 2. Comment on the relative purity of unknown A and how you determined this.
  - Unknown A is pretty pure. With a ramp rate of 1 °C/min, the change in temperature from when the sample started to melt to when it finished melting was only 0.9 °C.
- 3. What is the identity of unknown B and how did you determine this? Discuss both the IR and TLC data. Include a published IR spectrum in your report and comment on similarities between it and the spectrum you obtained.




Figure 2: trans-cinnamic acid IR spectra.

- Unknown B is trans-cinnamic acid. The experimentally obtained IR spectrum contains a distinct carboxylic acid peak (see Figure 2a), leading to the conclusion that B must be one of the seven possible compounds containing carboxylic acid functional groups. Comparison of these seven compounds with unknown B on the same TLC plates yielded identical and indistinguishable motion for unknown B and trans-cinnamic acid. Lastly, Figure 2 shows that the IR spectrum of unknown B and that of trans-cinnamic acid match nicely, especially in the fingerprint region.
- 4. Was the identity of unknown B definitive? In other words, were there any other possibilities?
  - Yes, it was definitive. Based on the exact match of the TLC dots (see notebook pages for sketch), unknown B is certainly *trans*-cinnamic acid.
- 5. What is the identity of unknown C and how did you determine this?
  - Unknown C is trans-stilbene. The IR spectrum obtained for the compound has a fingerprint region that matches nearly perfectly with the fingerprint region of the known IR spectrum of trans-stilbene.
- 6. Comment on the relative purity of unknown C and how you determined this.
  - Unknown C is quite pure. As per the above, the experimental and known IR spectra match nearly exactly.

Lab Assignment 1d CHEM 22000

## References

(1) SDBSWeb National Institute of Advanced Industrial Science and Technology https://sdbs.db.aist.go.jp/sdbs/cgi-bin/landingpage?sdbsno=1997 (accessed 10/20/2021).